


Lookup Tables

Image Processing & Analysis for Life Scientist

Olivier Burri, Romain Guiet & Arne Seitz

Look up tables

1. Create a new image
File→New→Image...
2. Recommended
Parameters:
Name: Ramp
Type: 8 bit
Fill with: Ramp
Width: 256 pixels
Height: 16 pixels
Slices: 1




Look up tables

3. Duplicate the image
Image→Duplicate...

OR

CTRL+SHIFT+D



4. Select the duplicated image
5. Invert the LUT
Image→Lookup
Tables
→Invert LUT
6. Align the images one beneath the other

Look up tables


Ramp (400%)

256x16 pixels; 8-bit; 4K



Ramp inverted (400%)

256x16 pixels; 8-bit (inverting LUT); 4K



7. Compare the pixel intensities at the same width position in the two images by moving the mouse pointer on the image.
8. Duplicate the first image again and explore another LUT.

Look up tables

Ramp (400%)

256x16 pixels; 8-bit; 4K



Ramp inverted (400%)

256x16 pixels; 8-bit (inverting LUT); 4K



Ramp Cyan (400%)

256x16 pixels; 8-bit; 4K



Ramp Fire (400%)

256x16 pixels; 8-bit; 4K

