


Working With Multichannel Images

Image Processing & Analysis for Life Scientist

Romain Guiet, Olivier Burri & Arne Seitz

Managing multiple channels on an image



1. **File** → **Open Samples** > HeLa Cells
2. “Split” the channels
Image → **Color** → **Split Channels**
3. Merge channels
Image → **Color** → **Merge Channels**
(Create Composite: OFF, Keep Source Images: ON)



Managing multiple channels on an image

1. **File** → **Open Samples** > HeLa Cells
2. “Split” the channels
Image → **Color** → **Split Channels**
3. Merge channels
Image → **Color** → **Merge Channels**
(Create Composite: OFF, Keep Source Images: ON)
4. **Image** → **Color** → **Merge Channels**
(Create Composite: **ON**, Keep Source Images: ON)

*Play with Brightness & Contrast on both images
(Try different color scenario).*



Managing multiple channels on an image

1. *File* → *Open Samples* > HeLa Cells
2. “Split” the channels
Image → *Color* → *Split Channels*
3. Merge channels
Image → *Color* → *Merge Channels*
(*Create Composite: OFF, Keep Source Images: ON*)
4. *Image* → *Color* → *Merge Channels*
(*Create Composite: ON, Keep Source Images: ON*)
5. **Select the RGB image**
Image → *Color* → ***Make composite (on RGB)***
Find the “Channels Tool...”
Set the mode to “**Composite**”, and change the display of channels.

