


Fourier Filtering

Image Processing & Analysis for Life Scientists

Olivier Burri, Romain Guiet & Arne Seitz


Getting the FFT

1. Open an image from HeLa_EdU
Duplicate the first channel only
2. Go to Process > FFT > FFT
3. Observe the FFT image




Modifying the FFT

1. Open an image from HeLa_EdU
Duplicate the first channel only
2. Go to **Process > FFT > FFT**
3. Observe the FFT image
- 4. Make a circle using
Process>FFT>Make Circular Selection**



Modifying the FFT


1. Open an image from HeLa_EdU
Duplicate the first channel only
2. Go to **Process > FFT > FFT**
3. Observe the FFT image
4. Make a circle using
Process>FFT>Make Circular Selection
5. Use **Edit>Selection>Clear Outside**
Make sure the background is black using the color picker




The Inverse FFT

1. Open an image from HeLa_EdU
Duplicate the first channel only
2. Go to **Process > FFT > FFT**
3. Observe the FFT image
4. Make a circle using
Process>FFT>Make Circular Selection
5. Use **Edit>Selection>Clear Outside**
Make sure the background is black using the color picker
- 6. Use Process>FFT>Inverse FFT**

We have applied a **Lowpass** filter: All high frequency components have been removed



Original




Modified FFT




Filtered

The Inverse FFT

1. Go back to the original image
2. **Process > FFT > FFT**
3. **Edit > Selection > Restore Selection**
4. Use **Edit > Clear**
5. Use **Process > FFT > Inverse FFT**



Original



Modified FFT





Filtered


We have applied a **Highpass** filter: All low frequency components have been removed

The FFT Bandpass Filter


1. Go back to the original image
2. **Process > FFT > Bandpass Filter**
3. Use settings as shown
4. Click OK
5. Use **Process>FFT>Inverse FFT**



We have applied a **Bandpass** filter: All periods between 5 and 10 pixels have been kept.



Resulting Filter



Resulting Image