



6.2.4. Regions of Interest - Practice

[Course](#) > [6. ROIs & Results](#) > [6.2. Regions of Interest \(ROIs\)](#) > Assessment

6.2.4. Regions of Interest - Practice Assessment

Region Of Interest

0 points possible (ungraded)

Which of the following statements are TRUE

- ☒ A Region Of Interest can be added by pressing the shortcut key "t" ✓
- ☒ You can define the fill color in the ROI properties ✓
- ☒ The button "Flatten" generates a RGB image containing with the ROI(s) ✓
- ☐ All the ROIs are saved automatically when an image is closed

Submit

i Answers are displayed within the problem

Operations on ROIs

0 points possible (ungraded)

Which of the following statements are TRUE

- ☐ the "OR" operator can be used to get the intersection of 2 ROIs
- ☒ the "AND" operator can be used to get the intersection of 2 ROIs ✓

☒ To measure multiples ROIs at once, you can Click "Deselect" then "Measure" ✓

☐ The function "Interpolate ROIs" does clever




Submit

i Answers are displayed within the problem

Of the use of ROIs

0 points possible (ungraded)

If needed, [download the HeLa cells images.](#)



and [download the macro "Detect_Cells.ijm"](#).

If dapi_chNbr = 1, please select recommended value(s) for measure_chNbr

☐ 1

☒ 2 ✓

☒ 3 ✓



If you set `measure_chNbr = 2`, and an several images from the the final "Results" table, the maximum value is 255, indicating that the image can be "saturated".

What would be your next decision :

☒ Discard this image from the dataset, and if too many images are saturated you discard the experiment ✓

☒ Look at the histogram of the channel to evaluate the "level" of saturation. ✓

☐ analyse individual cells and discard the cells containing "saturated pixels"



Considering the staining homogeneity in the different channels, you would rather:

☐ analyze the global values

☒ analyze individual cell values ✓



Submit

i Answers are displayed within the problem

[Learn About Verified Certificates](#)

© All Rights Reserved